

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 13.0cm (5.2")rotatable

for Mini vessel clips; TITAN



SKU: 244-680-13T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable