

YASARGIL APPLYING FORCEPS; RING HANDLE

rigid jaw; shaft 9.0cm (3.6")rotatable

for Standard vessel clips; TITAN



SKU: 244-640-09T
Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rigid , shaft rotatable