

YASARGIL APPLYING FORCEPS; RING HANDLE

jaw and shaft 11.0cm (4.4")rotatable
 for Standard vessel clips; TITAN



SKU: 244-642-11T
Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable , shaft rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable , shaft rotatable