

YASARGIL APPLYING FORCEPS; ROUND HANDLE

rotatable jaw; shaft 11.0cm (4.4")

for Standard vessel clips; TITAN

SKU: 244-637-11T
Categories: [Open Surgery](#), [Neurosurgery](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Applying Forceps for Aneurysm Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
SPECIFICATIONS	jaw rotatable