


YASARGIL MINI VESSEL CLIP

permanent; 7mm; opening 4.4mm

110 gms; curved



SKU: 244-467-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE722
AUTHORS	YASARGIL
Medicon	58.57.07
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE722
AUTHORS	YASARGIL
Medicon	58.57.07

SPECIFICATIONS

[curved](#), [permanent](#)