


YASARGIL MINI VESSEL CLIP

permanent; 4mm; opening 7mm

130 gms; angled



SKU: 244-476-04

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.43

SPECIFICATIONS [angled](#), [permanent](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.43

SPECIFICATIONS [angled](#), [permanent](#)