


YASARGIL MINI VESSEL CLIP

permanent; 7mm; opening 5.7mm

130 gms; angled



SKU: 244-476-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

Aesculap FE728

AUTHORS [YASARGIL](#)

Medicon 58.57.42

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE728

AUTHORS [YASARGIL](#)

Medicon 58.57.42

SPECIFICATIONS

[angled](#), [permanent](#)