


YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 4mm

110 gms; straight



SKU: 244-460-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE690

AUTHORS [YASARGIL](#)

Medicon 58.56.05

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE690

AUTHORS [YASARGIL](#)

Medicon 58.56.05

SPECIFICATIONS

[permanent](#), [straight](#)