


YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 4mm

110 gms; straight



SKU: 244-460-05

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE690
AUTHORS	YASARGIL
Medicon	58.56.05
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE690
AUTHORS	YASARGIL
Medicon	58.56.05

SPECIFICATIONS

[permanent](#), [straight](#)