


YASARGIL MINI VESSEL CLIP

temporary; 3mm; opening 3.3mm

90 gms; straight



SKU: 244-560-03

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE681

AUTHORS [YASARGIL](#)

Medicon 58.56.53

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE681

AUTHORS [YASARGIL](#)

Medicon 58.56.53

SPECIFICATIONS

[straight](#), [temporary](#)