


YASARGIL MINI VESSEL CLIP

temporary; 7mm; opening 4.6mm

70 gms; straight




SKU: 244-562-07

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE721
AUTHORS	YASARGIL
Medicon	58.56.58
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE721
AUTHORS	YASARGIL
Medicon	58.56.58

SPECIFICATIONS

[straight](#), [temporary](#)