


YASARGIL MINI VESSEL CLIP

temporary; 5mm; opening 3.8mm

70 gms; curved



SKU: 244-565-05

Categories: [Neurosurgery](#), [Open Surgery](#)


ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

SPECIFICATIONS [curved](#), [temporary](#)