


YASARGIL MINI VESSEL CLIP

temporary; 7mm; opening 4.4mm

70 gms; curved




SKU: 244-567-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE723
AUTHORS	YASARGIL
Medicon	58.57.58
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE723
AUTHORS	YASARGIL
Medicon	58.57.58

SPECIFICATIONS

[curved](#), [temporary](#)