


YASARGIL MINI VESSEL CLIP

temporary; 4mm; opening 3.5mm

80 gms; curved



SKU: 244-571-04

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.67

SPECIFICATIONS [curved](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.57.67

SPECIFICATIONS [curved](#), [temporary](#)