


## YASARGIL MINI VESSEL CLIP

permanent; 5mm; opening 3.8mm

110 gms; curved



**SKU:** 244-465-05

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE712
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.57.06
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>Aesculap</b>	FE712
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>Medicon</b>	58.57.06

## SPECIFICATIONS

[curved](#), [permanent](#)