


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 4mm; opening 3.6mm

110 gms; curved



SKU: 244-469-04T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT694T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.15
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT694T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.15

SPECIFICATIONS

[curved](#), [permanent](#)