


YASARGIL MINI VESSEL CLIP; TITANIUM

permanent; 7mm; opening 6mm

130 gms; straight



SKU: 244-475-06T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------------|--|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT717T |
| MATERIAL | Titanium |
| AUTHORS | YASARGIL |
| Medicon | 58.43.41 |
| SPECIFICATIONS | permanent , straight |



ZUSÄTZLICHE INFORMATIONEN

| | |
|-----------------|------------------------------|
| UNIT | pcs. |
| TYPE | Vessel Clips |
| Aesculap | FT717T |
| MATERIAL | Titanium |

AUTHORS

[YASARGIL](#)

Medicon

58.43.41

SPECIFICATIONS

[permanent](#), [straight](#)