


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 5mm; opening 4mm

70 gms; straight



SKU: 244-560-05T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT190T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.42.55
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT190T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.42.55

SPECIFICATIONS

[straight](#), [temporary](#)