


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 7mm; opening 4.6mm

70 gms; straight



SKU: 244-562-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT220T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.42.58

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT220T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.42.58

SPECIFICATIONS

[straight](#), [temporary](#)