


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 3mm; opening 3.2mm

90 gms; straight




SKU: 244-563-03T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.53
SPECIFICATIONS	straight , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.53

SPECIFICATIONS

[straight](#), [temporary](#)