


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 4mm; opening 3.6mm

80 gms; curved




SKU: 244-567-04T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.61
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.61

SPECIFICATIONS

[curved](#), [temporary](#)