


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 5mm; opening 4mm

70 gms; curved



SKU: 244-567-05T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.63
SPECIFICATIONS	curved , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.63

SPECIFICATIONS

[curved](#), [temporary](#)