


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 7mm; opening 6mm

70 gms; straight



SKU: 244-575-06T

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)


Aesculap FT217T

MATERIAL [Titanium](#)

AUTHORS [YASARGIL](#)

Medicon 58.43.70

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FT217T

MATERIAL [Titanium](#)

AUTHORS

[YASARGIL](#)

Medicon

58.43.70

SPECIFICATIONS

[straight](#), [temporary](#)