


YASARGIL MINI VESSEL CLIP; TITANIUM

temporary; 4.5mm; opening 6mm

70 gms; angled




SKU: 244-577-07T

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT306T
MATERIAL	Titanium
AUTHORS	YASARGIL
Medicon	58.43.73
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FT306T
MATERIAL	Titanium

AUTHORS

[YASARGIL](#)

Medicon

58.43.73

SPECIFICATIONS

[angled](#), [temporary](#)