


YASARGIL MINI VESSEL CLIP

permanent; 3mm; opening 3.3mm

110 gms; straight




SKU: 244-460-03

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE680
AUTHORS	YASARGIL
Medicon	58.56.03
SPECIFICATIONS	permanent , straight



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE680
AUTHORS	YASARGIL
Medicon	58.56.03

SPECIFICATIONS

[permanent](#), [straight](#)