


YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 8mm

110 gms; angled



SKU: 244-510-05

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	angled , temporary