


YASARGIL VESSEL CLIP; STANDARD

temporary; 17mm; opening 11.9mm

90 gms; angled



SKU: 244-510-17

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.17

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.59.17

SPECIFICATIONS [angled](#), [temporary](#)