


## YASARGIL VESSEL CLIP; STANDARD

permanent; 15mm; opening 8.7mm

200 gms; curved



**SKU:** 244-401-15

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>