


YASARGIL VESSEL CLIP; STANDARD

temporary; 10mm; opening 5.6mm

130 gms; angled



SKU: 244-513-10

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.10

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.10

SPECIFICATIONS [angled](#), [temporary](#)