


YASARGIL VESSEL CLIP; STANDARD

temporary; 11mm; opening 8mm

90 gms; straight



SKU: 244-521-11

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.41

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.63.41

SPECIFICATIONS [straight](#), [temporary](#)