


YASARGIL VESSEL CLIP; STANDARD

temporary; 12mm; opening 9.6mm

90 gms; straight; fenestrated Ø 3.5mm



SKU: 244-530-12

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.12

SPECIFICATIONS [straight](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.12

SPECIFICATIONS [straight](#), [temporary](#)