


YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 5.8mm

110 gms; angled; fenestrated Ø 3.5mm



SKU: 244-531-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.15

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.15

SPECIFICATIONS [angled](#), [temporary](#)