


## YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 6.5mm

110 gms; angled; fenestrated Ø 3.5mm



**SKU:** 244-531-07

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.17

**SPECIFICATIONS** [angled](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.17

**SPECIFICATIONS** [angled](#), [temporary](#)