


## YASARGIL VESSEL CLIP; STANDARD

temporary; 9mm; opening 9.1mm

100 gms; straight; fenestrated Ø 5mm



**SKU:** 244-536-09

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN


**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.39

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.58.39

**SPECIFICATIONS** [straight](#), [temporary](#)