


## YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 6.1mm

180 gms; curved



**SKU:** 244-403-09

**Categories:** [Neurosurgery](#), [Open Surgery](#)

### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>



### ZUSÄTZLICHE INFORMATIONEN

<b>UNIT</b>	<a href="#">pcs.</a>
<b>TYPE</b>	<a href="#">Vessel Clips</a>
<b>AUTHORS</b>	<a href="#">YASARGIL</a>
<b>SPECIFICATIONS</b>	<a href="#">curved</a> , <a href="#">permanent</a>