


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 6.1mm

180 gms; curved




SKU: 244-403-09

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	curved , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
SPECIFICATIONS	curved , permanent