


YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 7.9mm

90 gms; angled; fenestrated Ø 5mm



SKU: 244-540-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.65

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.65

SPECIFICATIONS [angled](#), [temporary](#)