


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 8.8mm

90 gms; angled; fenestrated Ø 5mm



SKU: 244-540-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.67

SPECIFICATIONS [angled](#), [temporary](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.58.67

SPECIFICATIONS [angled](#), [temporary](#)