


YASARGIL VESSEL CLIP; STANDARD

temporary; 7mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-541-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.07

SPECIFICATIONS [angled](#), [temporary](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.07

SPECIFICATIONS [angled](#), [temporary](#)