


YASARGIL VESSEL CLIP; STANDARD

temporary; 10mm; opening 5.8mm

110 gms; angled; fenestrated Ø 5mm



SKU: 244-543-10

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)


TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.20

SPECIFICATIONS [angled](#), [temporary](#)

Schließkraft
Closing Force



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.41.20

SPECIFICATIONS [angled](#), [temporary](#)