


## YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 7mm

110 gms; straight; fenestrated Ø 5mm



**SKU:** 244-545-05

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.25

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.25

**SPECIFICATIONS** [straight](#), [temporary](#)