


## YASARGIL VESSEL CLIP; STANDARD

temporary; 10mm; opening 8.8mm

110 gms; straight; fenestrated Ø 5mm



**SKU:** 244-545-10

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.30

**SPECIFICATIONS** [straight](#), [temporary](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**Medicon** 58.41.30

**SPECIFICATIONS** [straight](#), [temporary](#)