


YASARGIL VESSEL CLIP; STANDARD

temporary; 5mm; opening 5.6mm


110 gms; angled; fenestrated Ø 5mm



SKU: 244-551-05
Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.41.55
SPECIFICATIONS	angled , temporary



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.41.55
SPECIFICATIONS	angled , temporary