


## YASARGIL VESSEL CLIP; STANDARD

permanent; 15mm; opening 7.4mm

200 gms; curved



**SKU:** 244-403-15

**Categories:** [Neurosurgery](#), [Open Surgery](#)


### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**SPECIFICATIONS** [curved](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**AUTHORS** [YASARGIL](#)

**SPECIFICATIONS** [curved](#), [permanent](#)