


YASARGIL VESSEL CLIP; STANDARD

permanent; 7mm; opening 5.4mm

200 gms; curved



SKU: 244-405-07

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.55.37
SPECIFICATIONS	curved , permanent

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
AUTHORS	YASARGIL
Medicon	58.55.37
SPECIFICATIONS	curved , permanent