


YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 7mm

180 gms; straight



SKU: 244-400-09

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN


UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE750

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE750

AUTHORS [YASARGIL](#)

SPECIFICATIONS [permanent](#), [straight](#)