


## YASARGIL VESSEL CLIP; STANDARD

permanent; 9mm; opening 8.7mm

200 gms; angled



**SKU:** 244-410-09

**Categories:** [Open Surgery](#), [Neurosurgery](#)

### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE758

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.45

**SPECIFICATIONS** [angled](#), [permanent](#)



### ZUSÄTZLICHE INFORMATIONEN

**UNIT** [pcs.](#)

**TYPE** [Vessel Clips](#)

**Aesculap** FE758

**AUTHORS** [YASARGIL](#)

**Medicon** 58.55.45

## SPECIFICATIONS

[angled](#), [permanent](#)