


YASARGIL VESSEL CLIP; STANDARD

permanent; 20mm; opening 13mm

180 gms: angled




SKU: 244-410-20

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE798
AUTHORS	YASARGIL
Medicon	58.55.46
SPECIFICATIONS	angled , permanent



ZUSÄTZLICHE INFORMATIONEN

UNIT	pcs.
TYPE	Vessel Clips
Aesculap	FE798
AUTHORS	YASARGIL
Medicon	58.55.46

SPECIFICATIONS

[angled](#), [permanent](#)