


YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 6.2mm

200 gms; angled



SKU: 244-411-05

Categories: [Open Surgery](#), [Neurosurgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.50

SPECIFICATIONS [angled](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

AUTHORS [YASARGIL](#)

Medicon 58.55.50

SPECIFICATIONS [angled](#), [permanent](#)