


YASARGIL VESSEL CLIP; STANDARD

permanent; 5mm; opening 7.2mm

200 gms; curved



SKU: 244-417-05

Categories: [Neurosurgery](#), [Open Surgery](#)

ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE830

AUTHORS [YASARGIL](#)

Medicon 58.55.65

SPECIFICATIONS [curved](#), [permanent](#)



ZUSÄTZLICHE INFORMATIONEN

UNIT [pcs.](#)

TYPE [Vessel Clips](#)

Aesculap FE830

AUTHORS [YASARGIL](#)

Medicon 58.55.65

SPECIFICATIONS

[curved](#), [permanent](#)